



New Features

- **Mirroring in circumferential direction**
The parameter for „Mirroring in circumferential direction“ can now be set separately for each color. (E.g. usefull for counter printing.)
- **Storing of camera images**
New function for K500s which are equipped with CellGuard or HandyCam: The testcut cell images are stored now in the job ticket file. These can be displayed and printed by Helio-Disk. Thus an engraving log can be generated easily.
See illustration
- **Optimization of testcut layout and -data**
HelioKlischographs, equipped with CellGuard, perform an automatic testcut sequence. For setting up the testcut parameters several optimizations have been implemented. (see Read.me file)
- **Windows 2000 compatibility**
HelioDisk V. 2.33 ist prepared now for Windows 2000.

Hardware

No additional requirements.

Software

K500 version V. 1.24; K405 version V. 2.45